



Reference: PTXBC008764

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	kinesin-like protein KIF2C; CT139; KNSL6; MCAK; kinesin-like protein 6; mitotic centromere-associated kinesin; testis tissue sperm-binding protein Li 68n; kinesin family member 2C		

For research use only.

GENE DESCRIPTION	kinesin family member 2C
GENE ID	11004
ORIGIN SPECIES	Human
ACCESSIONS	BC008764
VECTOR	pDONR223
NCBI SYMBOL	KIF2C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]