



Reference: PTXBC008782

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NFKBIL2; IKBR; tonsoku-like protein; I-kappa-B-related protein; NF-kappa-B inhibitor-like protein 2; ikappaBR; inhibitor of kappa B-related protein; nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor-like 2; tonsoku-like, DNA rep		

For research use only.

GENE DESCRIPTION	nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor-like 2
GENE ID	4796
ORIGIN SPECIES	Human
ACCESSIONS	BC008782
VECTOR	pDONR223
NCBI SYMBOL	NFKBIL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr