



Reference: PTXBC008800

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDC2; CRCS10; MDPL; POLD; DNA polymerase delta catalytic subunit; CDC2 homolog; DNA polymerase subunit delta p125; polymerase (DNA directed), delta 1, catalytic subunit 125kDa; polymerase (DNA) delta 1, catalytic subunit; DNA polymerase delta 1, catalytic		

For research use only.

GENE DESCRIPTION	polymerase (DNA directed), delta 1, catalytic subunit 125kDa
GENE ID	5424
ORIGIN SPECIES	Human
ACCESSIONS	BC008800
VECTOR	pDONR223
NCBI SYMBOL	POLD1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]