



Reference: PTXBC008828

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZBTB22B; ZNF-X; zinc finger and BTB domain-containing protein 43; zinc finger and BTB domain-containing protein 22B; zinc finger protein 297B; zinc finger and BTB domain containing 43		

For research use only.

GENE DESCRIPTION	zinc finger and BTB domain containing 43
GENE ID	23099
ORIGIN SPECIES	Human
ACCESSIONS	BC008828
VECTOR	pDONR223
NCBI SYMBOL	ZBTB43
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]