

HMGA1-high mobility group AT-hook 1 Gene

Reference: PTXBC008832

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HMG-R; HMGA1A; HMGIY; high mobility group protein HMG-I/HMG-Y; high mobility group protein A1; high mobility group protein R; high-mobility group (nonhistone chromosomal) protein isoforms I and Y; nonhistone chromosomal high-mobility group protein HMG-I/H		

For research use only.

GENE DESCRIPTION	high mobility group AT-hook 1
GENE ID	3159
ORIGIN SPECIES	Human
ACCESSIONS	BC008832
VECTOR	pDONR223
NCBI SYMBOL	HMGA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtgcgagctcgaagtcagccagccctggcctcaagcaggaaaaaggacggcactgagaagcggggccggggcaggccgcgcaagcagcctccggtgagtcgggggacagcgtgtagggagtcagaaggagccagcgaagtccaacacctaagagacctc
ggggccgaccaagggaagcaaaaacagggtgctgccaagaccggaaaaccaccacaactccagggaaggaaccaaggggcagacccaaaaactggagaaggaggaagaggaggcatctcgaggagtcctcggaggaggagcagtga