

C1orf2-chromosome 1 open reading frame 2 Gene

Reference: PTXBC008854

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf2; COTE1; protein FAM189B; family with sequence similarity 189 member B		
For research use only.			

GENE DESCRIPTION	chromosome 1 open reading frame 2
GENE ID	10712
ORIGIN SPECIES	Human
ACCESSIONS	BC008854
VECTOR	pDONR223
NCBI SYMBOL	C1orf2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgccctgcctagtgaactcagccgctcgtgaccagccggccagcaccaggggcttaccacacccgcctccaccgacctgctgacggccctgcttacgctggggctgggtccaaagtgcctcctgggcatcctggtggtcaccctcagcatggtggcctcttcctgaccacc
accgagagcatcaagaggtcctgcccgtcttgggtgggttcgaacctgctcttcagcgtctgtgggtcaccattgtgccgctataatctgtacactctctgctattgctgctgcatccaaatctctccctggacctgctgcatagcgtggccctgagcgggtcagctcagggccactgg
gacctctgggctgcagctccccgccccagccccctctctacacacacatgctggacctggaggaattgtcccgcctgtgccccaccgccccctactatccccagagtatacctgcagctcagaacagatgcacagagcatcacgtacaatggctccatggacagccagctgccctgta
ccctaccgattgccccctcttatgaggcagtcagtgaggactacgaggaggaacaaattctgtccccctcagaatgagagtggtcttctgattgcaagggcactatggtcagggcaaggcatggcccaggtgtttaagtacaggggtgacgtgtgcctatgcaatgggggtgtaa