

SIX1-SIX homeobox 1 Gene

Reference: PTXBC008874

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	homeobox protein SIX1; BOS3; DFNA23; TIP39; sine oculis homeobox homolog 1; SIX homeobox 1		

For research use only.

GENE DESCRIPTION	SIX homeobox 1
GENE ID	6495
ORIGIN SPECIES	Human
ACCESSIONS	BC008874
VECTOR	pDONR223
NCBI SYMBOL	SIX1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgtatgctccgtctgtttgctttacgcaggagcaagtggcgtgctgtgctgaggttctgcagcaaggcggaaacctggagcgctgggcaggttctgtgtcactgcccgcctgcgaccacctgcacaagaacgagagcgtactcaaggccaaggcgggtgctgcctccacc
gcggcaacttccgtgagctctacaagatcctggagagccaccagttctgcctcacaaccaccccaactgcagcaactgtggtgaaggcgcattacgtggagccgagaagctgtgcggccgacctggggcgcgtgggcaaatatcgggtgcgcgaaaatttccactgcccgc
caccatctggacggcgaggagaccagctactgcttcaaggagaagtcgaggggtgtctctgctgggagtggtacgcgcacaatccctaccatcgccgctgagaagcgggagctggccgagggccaccggcctcaccaccacccaggtcagcaactggtttaagaaccggaggcaa
agagaccggcgccgaggccaaggaaagggaacaccgaaaacaataactctctccaacaagcagaaccaactctctctctggaagggggcaagccgctcatgtccagctcagaaggaattctcactccccaaagtccagaccagaactcggtcctctgctgcagggc
aatatgggccacgcccaggagctcaactattctctccgggcttaacagcctcgcagccagtcacggcctgcagaccaccagcatcagctccaagaactctctgctcgcccccctcactccagctctggtggacttgggtgctctaa