

EXOSC3-exosome component 3 Gene

Reference: PTXBC008880

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-102; PCH1B; RRP40; Rrp40p; bA3J10.7; hRrp-40; p10; exosome complex component RRP40; exosome complex exonuclease RRP40; ribosomal RNA-processing protein 40; exosome component 3		

For research use only.

GENE DESCRIPTION	exosome component 3
GENE ID	51010
ORIGIN SPECIES	Human
ACCESSIONS	BC008880
VECTOR	pDONR223
NCBI SYMBOL	EXOSC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgaacctgcgtctgcgcggctgaatctcgcgggcagcagggcgcgctgcacgcacagtactagtcaggtggtgtctccgggtgaggagctgctcctgccggaacaggaggacgcggaaggccctgggggtcagtgagcgaccgttgagcctgaatgctagag
cgtgctcgcgggtgcgcgtgtatgcggtccggccttcggcgtgtggggaccgctgctgtgacccaagtgcggccgctccgtcacaaggagcccgcagtggcagcggcgccgggtttactgggtgactctcagcagaagcggtatgtccagttaaaaggagaccatgtgat
tggcatagtgcagctaaatctggagatatacaaaagtgtggtggaggagtgagccagctcttctgtctactgtcattgaaggtgcaactaaaagaaacagacaaatgtgcaggttgagatctcatatggccagtttgggtgtaataaagacatggaaccagagatggtctgt
attgacagctgtggacgagccaatggaaatgggtgctcattggacagatggtctgttttaaaagtgaactctgggttaattagaagctattagctccagattgtgaaatcatacagggaagtgggaaactctatccactggagatagtattggaatgaatggaagaatatgggttaaggcaaa
aaccatccagcagactttaatttggcaacatttagaagcttgaacacatgacgtcagatcaaaagaaacagatctctccagattggcagaagtga