

MID1IP1-MID1 interacting protein 1 (gastrulation specific G12 homolog (zebrafish)) Gene

Reference: PTXBC008908

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	G12-like; MIG12; S14R; STRAIT11499; THRSPL; mid1-interacting protein 1; MID1 interacting G12-like protein; MID1 interacting protein 1 (gastrulation specific G12-like); gastrulation specific G12 homolog; spot 14-R; spot 14-related protein; MID1 interacting		

For research use only.

GENE DESCRIPTION	MID1 interacting protein 1 (gastrulation specific G12 homolog (zebrafish))
GENE ID	58526
ORIGIN SPECIES	Human
ACCESSIONS	BC008908
VECTOR	pDONR223
NCBI SYMBOL	MID1IP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgcaaatctgcgacacacacaccagaagcactcgtctttaacgcatgaatcgttcattggcgccgtgaacaacatggaccagacggtgatgggtgccagcttgcgcgacgtgccccggctgacccgggttagacaacgatgttggtggaggtaggcggcagtggc
ggctgcctggagagcgcacgccccagtcctcggggaagcgcaatggcagcttttcgcccctcgggacatgtacagccactcgtcttcaagtcacccgcaacgacatcagtggtgggtcctgcaccagcgcctccaccggctgggagcgaggaggcagtg
gcctggaagtccaaggacatcctggtggacctgggcccattggagggtgcggacgcccgcgaagaagacctggaacagcagttccactaccactgcgcgggtgcacactgtgctctcgaaactacgcgcaaaagccaacatcctactaacagatacaagcaggagatcggttc
ggcaattggggccactga