



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
 contact@proteogenix.fr
 www.ProteoGenix.science

FAM3A-family with sequence similarity 3, member A Gene

Reference: PTXBC008912

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM3A; DLD; DXS560S; XAP-7; cytokine-like protein 2-19; family with sequence similarity 3 member A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 3, member A
GENE ID	60343
ORIGIN SPECIES	Human
ACCESSIONS	BC008912
VECTOR	pDONR223
NCBI SYMBOL	FAM3A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggtgtggcagcccccctgcgcattggtgctcagtcgtcattggtgtggtgtgcacatggtatcgtggtcacatgcctctctgggtggcgtggcagatggctttctccgcatccagcaactcttcaccagtcacagagagctcgtgactcagcggccacggggcagggaagtaacagtggtgctt
gccccagcgtctgctcgtgaggagcacctggcctccgcgtggtgcagcggggccgcacacgtcattgggtcccaagatctgctgtaggacagatgctgattgagcagcgtcagaagcaacatggggccgcgggctgaacatcgccctggtgaacggggctcagcggcgagctcattcag
gcggcggcctttgacatgttggcggcgagatgcacagcagctgattgatttatttgcggccatcgacagcagaagcctgctgttcgtgcctatctacagcaccagccacccaagatgaatgaagacacagagaagctcttcagctgagctcgtggcagcagacccaagggagcctggcctt
ccggcagcgtgtgtgttttctcggcccaagcgtgtgcagacaagaagcccccttcagcaacatctgaagaacatgaagcagcagatgaagcctgcggccagcgtcctcgaagatgaagcgtctatccgcgaaagcaacacacacagcag

More info: Please contact us for any question: contact@proteogenix.fr