

TXNDC12-thioredoxin domain containing 12 (endoplasmic reticulum) Gene

Reference: PTXBC008913

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AG1; AGR1; ERP16; ERP18; ERP19; PDIA16; TLP19; hAG-1; hTLP19; thioredoxin domain-containing protein 12; ER protein 18; ER protein 19; anterior gradient homolog 1; endoplasmic reticulum protein ERp19; endoplasmic reticulum resident protein 18; endoplasmic		

For research use only.

GENE DESCRIPTION	thioredoxin domain containing 12 (endoplasmic reticulum)
GENE ID	51060
ORIGIN SPECIES	Human
ACCESSIONS	BC008913
VECTOR	pDONR223
NCBI SYMBOL	TXNDC12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagacgcccgcctcgtctcggggccacctgttgcctggcctcagttctcctcgtcatctctctgatggacataatggccttgaaagggtttggagatcatattcattggagacactggaagatgggaagaagaagcagctgccagtggactgcccctgatggtgattatca
taaatcctggtgtggagcttgcaaaagctctaaagcccaattgcagaatctacggaaattcagaactctccataattttgtatggtaaatcttgaggatgaagaggaaacccaacatgaagatttcagccctgacgggggttatattccacgaatcttttctggatccagtggaagggtg
catcctgaaatcatcaatgagaatggaaacccagctacaagtattttatgtcagtgccgagcaagttgtcaggggatgaaggagctcaggaaggctgacgggtgatgccttcagaagaacatcttgaagatgaattgtaa