



Reference: PTXBC008957

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C19orf6; ASBABP1; MBRL; R32184_3; aspecific BCL2 ARE-binding protein 1; transmembrane protein 259		

For research use only.

GENE DESCRIPTION	chromosome 19 open reading frame 6
GENE ID	91304
ORIGIN SPECIES	Human
ACCESSIONS	BC008957
VECTOR	pDONR223
NCBI SYMBOL	C19orf6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]