

ATP6V1E2-ATPase, H⁺ transporting, lysosomal 31kDa, V1 subunit E2 Gene

Reference: PTXBC008981

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP6E1; ATP6EL2; ATP6V1EL2; VMA4; V-type proton ATPase subunit E 2; ATPase, H ⁺ transporting, lysosomal 31kDa, V1 subunit E2; V-ATPase subunit E 2; testis secretory sperm-binding protein Li 235P; vacuolar proton pump subunit E 2; vacuolar-type proton-trans		

For research use only.

GENE DESCRIPTION	ATPase, H ⁺ transporting, lysosomal 31kDa, V1 subunit E2
GENE ID	90423
ORIGIN SPECIES	Human
ACCESSIONS	BC008981
VECTOR	pDONR223
NCBI SYMBOL	ATP6V1E2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccctgagtgatgctgatgtgaaaaagcagattaagcaccatgatggcttcattgagcagggaagccaatgagaaagcagaggaaatcagccaaggctgagggaaggttaacattgagaaaggacgcctcgtgcaaaccaacgactgaagattatggagtattatgagaaaaag
gagaagcagatagagcagcagaagaataatcctgatgtccaccatgagggaatcaggcgaggctgaaagtctgagagccggaatgacctatctcagatttgcctcagtgaggcgaagctgagactcagcaggatttgaggagaccagaggctaccaggggctgctgataaactg
gtgtccagggtctgtccgactgctggaacctgtgatgtgtacgtgcccggccacaagacctctcctgtgtgaggctgtgtacaaaaagccatccccgagtacatgacaatttccagaaacatgtggaggtccagattgataaagaggcaccctggctgtgaatgcagctggag
gtgtggaggtctacagtggaatcagagaataaaggttcaataaccttggaaagccgactggatctctcagccaagcaaaagatgccagaatacgaatggccttgtttggtgtaacaccaacagaagttctttatataa