



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C1QC-complement component 1, q subcomponent, C chain Gene

Reference: PTXBC009016

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1Q-C; C1QG; complement C1q subcomponent subunit C; complement component 1, q subcomponent, C chain; complement component 1, q subcomponent, gamma polypeptide; complement C1q C chain		

For research use only.

GENE DESCRIPTION	complement component 1, q subcomponent, C chain
GENE ID	714
ORIGIN SPECIES	Human
ACCESSIONS	BC009016
VECTOR	pDONR223
NCBI SYMBOL	C1QC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtacgtggggccagctccctgccacactgggtcgaagctgctgctctgctgctgctccctcagagggccgaagccaacacagcgtctacgggatccagggaftccggcctccttggggcaccagggaagatgggtgtacgacggactcggggggcccaaggggg
agccaggaaatccagccattccgggatccgaggaccacaaggcgagaagggaagaaccggcttaccggccatctcgggaaaaatggcccattggacccttcgggatccagggggtgccggcccatgggcatccctggagagccaggtgagggggcatagataaagcag
aaattccagtcagtgttaccggctactcggcagaccacagccccctgaccacaacagcctgatcagattaacgcgttctcacaaccgcagggagatattgacacagacgtctgcaagtactctgcaaatcccgccctctactcttctgtacacacgcgtgcatacagcca
acctgtgcgtctgctgactccagcagcgtcgaagtggatcctctgtggccacatgctcaaaacaaatcagctcaactggcggtgtgctgctgagggtgcaggtggggcgaggaggtgtggtgctgctgctcaatgactactacgacatgtggggcatcagggtcgtgacagcgtct
ctccgcttctcgttccccaagcag

More info: Please contact us for any question: contact@proteogenix.fr