

PPP2R3B-protein phosphatase 2 (formerly 2A), regulatory subunit B'', beta Gene

Reference: PTXBC009032

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NYREN8; PPP2R3L; PPP2R3LY; PR48; serine/threonine-protein phosphatase 2A regulatory subunit B'' subunit beta; NY-REN-8 antigen; PP2A, subunit B, PR48 isoform; protein phosphatase 2 (formerly 2A), regulatory subunit B'', beta; protein phosphatase 2, regula		

For research use only.

GENE DESCRIPTION	protein phosphatase 2 (formerly 2A), regulatory subunit B'', beta
GENE ID	28227
ORIGIN SPECIES	Human
ACCESSIONS	BC009032
VECTOR	pDONR223
NCBI SYMBOL	PPP2R3B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacctggacggggacggcgccctgtccatgttcgagctcgagtcgacttctacgaggagcagtgccgaaggctggacagcatggccatcgaggccctgcccctccaggactgcctctgccagatgctggacctgggtcaagccgaggactgaagggaagatcacgctgcaggacctg
aagcgctgcaagctggccaacgtcttcttcgacaccttctcaacatcgagaagtacctcgaccacgagcagaaagagcagatctccctgctcaggagacgggtgacagcgcgccgccccgagctctcgactgggagaagtagcgggccgaggagtagacatccctgtggccgagga
gaccgcgggagagccctgggaggacgggttcgaggccgagctcagccccgtggagcagaagctgagtgcgctcgctccccgctggccagaggcccttcttcgaggcgccctaccgctggcgccgtggacctgtacgagtacgctgccccgaggacgagacctggagccgc
tgtga