

LIN37-lin-37 homolog (C. elegans) Gene

Reference: PTXBC009071

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	F25965; ZK418.4; lin-37; protein lin-37 homolog; antolefinin; lin-37 homolog; protein F25965; lin-37 DREAM MuvB core complex component		

For research use only.

GENE DESCRIPTION	lin-37 homolog (C. elegans)
GENE ID	55957
ORIGIN SPECIES	Human
ACCESSIONS	BC009071
VECTOR	pDONR223
NCBI SYMBOL	LIN37
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcctgtgaagtgaaagtggagaatcagagctggagatggccaaagccgggaaccaactggatgtgtcttcagctgtctgtggagaaggtcacatggacaggagcgtctggatgaggaagctgggaaaacaccctcagacaccacaataaggactgtccatgcag
ccactggcaaaaggccatctgcccgttccccaccagcggaggaagaagaggaggatggatgagggctgctgaggaggccgcagcgaaccaacacatatgtgatcaagctgttcgaccggagcgtggacttggccagttcagcgagaacacgccactgtacccaatc
tgccgcctggtgagcaacagccctctgtgctgagcgtgaatgtctccagctcaccctgccccgctgctgaggatgaggaggctcagaggttaaccaacagcaagagtcgtgatgtgtacaagctgcccaccacacccccggggccaccggagatgcctgcag
atccgcataccatctccactgcagcctgagatgcagggcaccctgacgatgagccctctgagcccgagccctcaccctccacactcatctatcgcaacatgcagcgtggaaacgcatcccgagaggtggaaggaggcctctatcggaaccagcttcgttactcagaagcatga
agatcctacgagatgtacgaacgacagtga