

## FAM181A-family with sequence similarity 181, member A Gene

Reference: PTXBC009073

|                               |                                                                          |                     |                |
|-------------------------------|--------------------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                                      | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                                     | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | protein FAM181A; C14orf152; family with sequence similarity 181 member A |                     |                |
| <b>For research use only.</b> |                                                                          |                     |                |

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | family with sequence similarity 181, member A                                                     |
| <b>GENE ID</b>          | 90050                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC009073                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | FAM181A                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggcatccgacagtgatggaagatgctgctgaactctggaactctggcgtccagcgacatcaaggcagccctggataagtcgcacccctgccgccgctccgtggaccatcgcaagtacctgcagaagcagctcaagcgtctccagaagtattccggctcccggggccttct  
ggcagagctgctgagccctacctgaaaaaggggctgaggaccggccaggagctgctctggaattggccctgattccagcccgccgggggtggggctgcaaggagaaggtgctgaggaaacctacaggaggaatgtctgtaaggagcagctcccacagaggcagc  
atccagaagctgccagcctggccaggtgcccagtgaggaaaagacagctgcccgttctctgggaagagccaaggccacccacagctaccatgtgggctggagggggactggggcccgaggaggacccctctatgagggttaagaaaaattgcaagggttgaggccct  
gggacctgagactacccctgtgtccatgtctccaaggccctggctgaaaaggagccgctcaagatgctgggtctcttgggtggccgcgtcaatgctgaggtgctgccccttcagtaaccatggacagcccatctatccggggcccttggggcactgctcagagtcctgtccc  
cagcctgggcttgggaggaagagcccagccttcccggggagctggcgacactctgcaaggatgtggacggcctggggcagaaaggtgtgagcccggtgtgtgctgaaccatcccaccaagcagcgtgccccaccatcttcaatgtcttggctacctctag