

FAM73B-family with sequence similarity 73, member B Gene

Reference: PTXBC009114

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM73B; FAM73B; C9orf54; mitoguardin 2; family with sequence similarity 73, member B		

For research use only.

GENE DESCRIPTION	family with sequence similarity 73, member B
GENE ID	84895
ORIGIN SPECIES	Human
ACCESSIONS	BC009114
VECTOR	pDONR223
NCBI SYMBOL	FAM73B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgacagcctgatgaccaaggctgagaagagccccaaaggcttctggagagctacgaggagatgctgagctatgccctgcccagacctggccacaacacacggctggagctggagggccgaggggtggtatgcatgagcttctgacatcgtgctggacttcctca
tggagccttcgaggacctggagaacctccggcctcggctgctgcccgtcctgctggaaccgctggctgtcagacagctcaaggagacggccttggccactgcttggctgctgctgaaagccaagaggagggctgctgctgctgctgcttctactcctgacg
gagcatgtgagccctgctagccttcggctcctggacccaagcctcagcttgcgaagtctgtcttctcaagcaccagatgtgctgctgctgagggacatgttgacctggacaatgtgctgctacacgtcactgcccgcgctggcagacacatcctgagctgtccggcgccgc
agcgagatattgctgggtacctgggggtgcccgcgcccagcagcgcagggcgtgaatggggcgtgccccgagagaatggcccctgggggagctgcagtag