

TSPYL4-TSPY-like 4 Gene

Reference: PTXBC009116

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	dJ486I3.2; testis-specific Y-encoded-like protein 4; TSPY-like protein 4; TSPY like 4		
For research use only.			

GENE DESCRIPTION	TSPY-like 4
GENE ID	23270
ORIGIN SPECIES	Human
ACCESSIONS	BC009116
VECTOR	pDONR223
NCBI SYMBOL	TSPYL4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggactcactggagccatcatcaagagttgtcaaacgtaaatgccaggctgacagggccttctcagcttgagcgcaagttggccgcatgcgaagctccacatgcagcgcagaagtttcattatccagaatatccagggtttctgggttactgccttgcgaaccacccccagctgt
cacctatgatcagtgcccaagatgaagacatgctgaggtacatcatcaatgtggaggtggagggagcttaaacccccagcagagctgcaaatcaagttcatctttcaggcgcaaccctacttccgaatgaggggctgtgcaagggaatgaacgcagatcctctggccgggtgtgtct
ctttccactccaatccgtggcaccgagggccaagacccccaggtcatatccacagaaacgggaagggaacactatccctagtttctcaactggttttcagaccacagccttctagaattgcagagaattgcagagattatcaaggagaactgtggcccaatcccctacaatactactga
tgggtgaaggggcccgtagaggaaatcgagggcccaaggcagccagtgagagcgccagatccttcaggttcagttcgctaa