

C14orf159-chromosome 14 open reading frame 159 Gene

Reference: PTXBC009182

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0317 protein C14orf159, mitochondrial; chromosome 14 open reading frame 159		
For research use only.			

GENE DESCRIPTION	chromosome 14 open reading frame 159
GENE ID	80017
ORIGIN SPECIES	Human
ACCESSIONS	BC009182
VECTOR	pDONR223
NCBI SYMBOL	C14orf159
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccttcacactccacgtgaggtccgccctccctctgccaataaggagtttgattctacaaaagaaccaacatcagaataacatccagcatgggtggagagctccgaccagccagcctgggtgctctgccaggtcccttgctccagctttgaaagatttgcaggtcaacactggc
ctctacccctgctggccagagtgagccagaaaagtgatgctgccccctcaagggtctatctcagagaccagcctcaggaccttcacagtgctggcagagaacctcactcatgactgttgcaggaacgagctgccccctccactgttctgactcttccctatcacctggcgcctgt
cctggctgagatcttggatcatggagaagtcagccctggagaatccactggctctgccatctgtccactccatcgcccaagccagaacccgcctgtacaggattaa