

CLTA-clathrin, light chain (Lca) Gene

Reference: PTXBC009201

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LCA; clathrin light chain A; clathrin, light polypeptide (Lca)		
For research use only.			

GENE DESCRIPTION	clathrin, light chain (Lca)
GENE ID	1211
ORIGIN SPECIES	Human
ACCESSIONS	BC009201
VECTOR	pDONR223
NCBI SYMBOL	CLTA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgagctggtacgcttgcggcggcctgcccggcgccctggcggtcccgctggggaacggagtgccggcgccggcgaagaagaccggctgcggccttcttggcgagcaagagagcgagattgcggcatcgagaacgacgaggccttgcctcctgacggcgg
cgccccgggccccagccgacggcgagccgggggggtccggatgctgtgagtgatgaatgggtgaatactaccaggaagtaattgccaacagacagttatgcagctatttcacaagtgatcgattgcagtcagagcctgaaagtatccgtaaatggagagaagaacaa
atggaacgcttggaagcccttgatgccaattctcggaagcaagaagcagagtggaagaaaaaggcaataaaggagctagaagaatggtatgcaagacaggacgagcagctacagaaaacaaaagcaacaacaggcgagcagaagaagcctttgtaaatgacattgacgagtcgic
cccaggcactgagtggaacgggtggccggctgtgtgactttaacccaagctagcaagcaggccaagaatgtctcccgatgcgtcagtcctctcctcaagcaggccccgctgtgcactga