

## PAGE5-P antigen family, member 5 (prostate associated) Gene

Reference: PTXBC009230

|                |                                                                                                                                                                                                                                                                |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                            | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                           | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CT16; CT16.1; CT16.2; GAGEE1; PAGE-5; P antigen family member 5; P antigen family, member 5 (prostate associated); cancer/testis antigen 16.1; cancer/testis antigen family 16, member 1; cancer/testis antigen family 16, member 2; g antigen family E member |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | P antigen family, member 5 (prostate associated)                                                  |
| <b>GENE ID</b>          | 90737                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC009230                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | PAGE5                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgcaggcgccatgggcccgttaacgtggctgggctgggaacgaggaggaagtgagagatatgagtgagcatgaacaagatcccaatctcagaaaggagaaatgaccaagagctctccagccagttggacctgtgattgtccagcagccactgaggaaaaacgtcaagaagaggaaccaccaactgataatcagggtattgcacctagtgaggagatcaaaaatgaaggagcacctgctgtcaagggaactgatgtggaagctttcaacagggaactggctctgcttaagatagaggatgcacctggagatggctctgagtcaggaggaggactctgccactttgatccacataaagtctggaagcaggtgaagggaactatag