



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

KPTN-kaptin (actin binding protein) Gene

Reference: PTXBC009249

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	2E4; MRT41; kaptin; actin-associated protein 2E4; kaptin, actin binding protein		
For research use only.			

GENE DESCRIPTION	kaptin (actin binding protein)
GENE ID	11133
ORIGIN SPECIES	Human
ACCESSIONS	BC009249
VECTOR	pDONR223
NCBI SYMBOL	KPTN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr