



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TMEM203-transmembrane protein 203 Gene

Reference: PTXBC009283

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

HBEBP1; transmembrane protein 203; HBeAg-binding protein 1

For research use only.

GENE DESCRIPTION	transmembrane protein 203
GENE ID	94107
ORIGIN SPECIES	Human
ACCESSIONS	BC009283
VECTOR	pDONR223
NCBI SYMBOL	TMEM203
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctctctcgtcccccggagctcgtgacgctgctgcacactcttcgagacatcttcgtcaactctgtggccctctgttggtctctctgtgctcgtggccactcgtgtggatggccgtgccccggggcctctcctgttggaacgtgttctgtgctttcttcgcgcgtgaacgggctcagcaactac
ttcacaccacatctgtccctgttcctcttcaggagtcagagagaagcggctgcgcggtgtccgcctttcttggtgactctacggctcctgagttcaagttcgtctcagatgctgttggtccagaagctggcgagagacactcgggagctctggttcggcctcattacgtccccgcttctatctc