

GABARAPL1-GABA(A) receptor-associated protein like 1 Gene

Reference: PTXBC009309

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GABARAPL1-a; APG8-LIKE; APG8L; ATG8; ATG8B; ATG8L; GEC1; gamma-aminobutyric acid receptor-associated protein-like 1; GABA(A) receptor-associated protein-like 1; GEC-1; early estrogen-regulated protein; glandular epithelial cell protein 1; GABA type A rece		

For research use only.

GENE DESCRIPTION	GABA(A) receptor-associated protein like 1
GENE ID	23710
ORIGIN SPECIES	Human
ACCESSIONS	BC009309
VECTOR	pDONR223
NCBI SYMBOL	GABARAPL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagttccagtacaaggaggaccatccctttgagatcggaaaaaggaaggagaaagatccggaagaaatccggacaggggtcccggtgattgtagagaaggtccaaaagccaggggtgcctgatctggacaagaggaaagtagtgcctctgacctactgttggccagttc
tacttcttaalccggaagagaatccacctgagacctgaggagcgttattcttcttgcacaacaccatccctcccaccagtgtctacatgggcccaactgtatgaggacaatcatgaggagactatttctgtatgtggcctacagtgatgagagtgtctatgggaatga