

PDRG1-p53 and DNA damage regulated 1 Gene

Reference: PTXBC009334

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf126; p53 and DNA damage-regulated protein 1; p53 and DNA damage regulated 1		

For research use only.

GENE DESCRIPTION	p53 and DNA damage regulated 1
GENE ID	81572
ORIGIN SPECIES	Human
ACCESSIONS	BC009334
VECTOR	pDONR223
NCBI SYMBOL	PDRG1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctatcacccgagcagagcagtgctgcgtacctgtagaagtgaggagctgccgaggaggtgctggcggaacagcgcagattgtggacctggacactaaaaggaatcagaatcagagggcctgagggccctgcagaaggatctcagcctctctgaagatgtgatggtt
gcttcgggaacatgtttatcaagatgcctcacctgagacaaaggaaatgattgaaaaagatcaagatcatctggataaagaataagaaaaactgcggaagcaacttaagtgagggtcaaccgctttttgaggcccaaggcaaacggagctgaagggttttaactgaacccctcaac
caggatgagcttaagctctcaaggctcatctgaaaggatga