

RGS10-regulator of G-protein signaling 10 Gene

Reference: PTXBC009361

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	regulator of G-protein signaling 10; regulator of G-protein signalling 10		
For research use only.			

GENE DESCRIPTION	regulator of G-protein signaling 10
GENE ID	6001
ORIGIN SPECIES	Human
ACCESSIONS	BC009361
VECTOR	pDONR223
NCBI SYMBOL	RGS10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcaaccgcgcgtgagccggctgagcaggaagcggccgctcagacatccacgacagcgatggcagttccagcagcagccagagcctcaagagcacagcacaatggcgccatccctggagaatctgctggaagaccagaaggcgtgaaagatttagggaatttt
aaaaaagggaattcagtgaaagaaatgtttgtttggctagcatgtgaagattttaagaaatgcaagataagacgcagatgcaggaagggcaagagagatctacatgacctttctgtccagcaaggcctcatcacaggtaacgtggagggcagctcggctcaacgagaagatcctgg
aagaaccgcaccctctgatgttccagaaactccaggaccagatcttaatctcatgaagtacgacagctacagccgcttcttaagctctgactgtttttaaacaacaagcgaaccgaggaagaggaagaagatttcctgatgctcaactgcagctaaaagagcttccagaattataacaca
tga