

STMN3-stathmin-like 3 Gene

Reference: PTXBC009381

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SCLIP; stathmin-3; SCG10-like protein; stathmin-like 3; stathmin 3		
For research use only.			

GENE DESCRIPTION	stathmin-like 3
GENE ID	50861
ORIGIN SPECIES	Human
ACCESSIONS	BC009381
VECTOR	pDONR223
NCBI SYMBOL	STMN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccagcaccattccgctacaaggagaagatgaaggagctgtcgggtgctgtcgtcatctgtcctctcttctacacacagccgcaccccaataccgtctaccagtacggggacatggaggtgaagcagctggacaagcgggcctcaggccagagcttcgaggtcctcctcaagtc
cccttctgacctgtcccgagagccctatgtctcctccccacccaagaagaaggacacctccctggaggagctgcaaaagcggctggaggcagccgagggaggaagacgcaggagcgaggtgctgaagcagctggcgagcgccgcgagcagcgcgaggtg
ctgcacaaggcgtggaggagaatacaactcagcccgaggcggaggagaagctcaactacaagatggagctcagcaaggagatccgcgagggcacacctggccgactgcgcgagcggtgctgcgagaaggagctgcacgcggccgaggtgctgcaggagaagagcag
cgagaagagatgtcgggctaa