

FAM163A-family with sequence similarity 163, member A Gene

Reference: PTXBC009382

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM163A; A230106N23Rik; cebelin; family with sequence similarity 163, member A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 163, member A
GENE ID	148753
ORIGIN SPECIES	Human
ACCESSIONS	BC009382
VECTOR	pDONR223
NCBI SYMBOL	FAM163A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacagcgggaacgggtgtgatcactggcgggaatcctagctacgggtgacctcctctgcattgccgtcgtgtgactgcaggctccagtattactgctgcaagaagagcggaaaccgaggttgacagcaggaggaggaggcgggagcagacctcccacgcatccagaggccc
cactgcaatgcctgcagctcccaagccctggacggcagaggcagcctggcgctctcaccagcagccctgcagccagccctgtgggtggcgcgagccactgcactacgtctccccatacagctccccctttacatacggagcgtgacatggtgccaatgggggtggagg
cgagaggctctccttgcctccacatactacaagaggggggacccccatccctcaaatggcagcaccagagttaccgggtgacctggccaggctctggcgctgaggccttccaatccaagggtattagtagcagacgtgtaa