



Reference: PTXBC009387

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGAT; VAT1; VMAT1; chromaffin granule amine transporter; solute carrier family 18 (vesicular monoamine transporter), member 1; solute carrier family 18 (vesicular monoamine), member 1; solute carrier family 18 member 1; vesicular amine transporter 1; solu		

For research use only.

GENE DESCRIPTION	solute carrier family 18 (vesicular monoamine), member 1
GENE ID	6570
ORIGIN SPECIES	Human
ACCESSIONS	BC009387
VECTOR	pDONR223
NCBI SYMBOL	SLC18A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]