



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

AGFG2-ArfGAP with FG repeats 2 Gene

Reference: PTXBC009393

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

HRBL; RABR; arf-GAP domain and FG repeat-containing protein 2; HIV-1 Rev-binding protein-like protein; Rev/Rex activation domain binding protein-related; arf-GAP domain and FG repeats-containing protein 2; nucleoporin; rev/Rex activation domain-binding pr

For research use only.

GENE DESCRIPTION	ArfGAP with FG repeats 2
GENE ID	3268
ORIGIN SPECIES	Human
ACCESSIONS	BC009393
VECTOR	pDONR223
NCBI SYMBOL	AGFG2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggatgacggcgcgaagaaggcccccggccggcgccggcgccgacggcgccgaaggcggagcggcgccgcctcggaaggtgtgtgtgccgtcgggtgcggagctgggtgctgcacccagccgggaaccgcactgtctcagatgcgccagcgcggggtcaccta
cctggatataccgggtgcggacgttcgtgtgcacactctgcctgcgcctctgagaggctgacacccccctactcgtgtaacgtcaatctcaacaaatctcagcctgagagtagtatatactcaatccggaggaatgaggtgttcgccgaagatgttggtgggtctgttgatcgtcga
catcttgcacacagatgccagatccctcagaagaatgaagagatttcccaaaaatagagaagaagaatgccagacacgtcccaaggaggtgttcgccacattg

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1