

IL32-interleukin 32 Gene

Reference: PTXBC009401

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL-32alpha; IL-32beta; IL-32delta; IL-32gamma; NK4; TAIF; TAIFa; TAIFb; TAIFc; TAIFd; interleukin-32; interleukin-32 eta; interleukin-32 small; interleukin-32 theta; natural killer cell transcript 4; natural killer cells protein 4; tumor necrosis factor a		

For research use only.

GENE DESCRIPTION	interleukin 32
GENE ID	9235
ORIGIN SPECIES	Human
ACCESSIONS	BC009401
VECTOR	pDONR223
NCBI SYMBOL	IL32
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtctccgaaggtctctctgatgacatgaagaagctgaagggccgaatggtaatgctcctcctactctgtctcaggggttggggcctgggtctcagctgtgacactgaggacactgtgggacacctgggacccctggaggacaaggatcggcccttgggtgccaactctgcc
tctcttcacagcaccagggccatagaagaatttatgataaaatgcaaaatgcagaatcaggagctggacaggtgatgtcggagcctggcagagctggaggacgacttcaagagggtacctggagacagtgccggcttattatgaggagcagcaccagagctcactctctactgaaaa
agaaagagatggattacgggtgccgagggcaacagatccctgtccggatgttgaggatcccgaaccgaggagcctggggagagctttgtgacaaggtcatgatgttccaggccatgctgcagcggctgcagacctgggtggcaggggttctggcctgggtgaaggagaaggt
gggtggcctggtccatgcagtcaggccctctggaacagttccagagttctgtctctgtcagagcttcatgtcctcttccagtcctacggagcccccacgggggacaaggagagctgacccccagaagtgtctgaacccaatctcaaaatga