

KIF26A-kinesin family member 26A Gene

Reference: PTXBC009415

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KIF26A variant protein; kinesin-like protein KIF26A; kinesin family member 26A		
For research use only.			

GENE DESCRIPTION	kinesin family member 26A
GENE ID	26153
ORIGIN SPECIES	Human
ACCESSIONS	BC009415
VECTOR	pDONR223
NCBI SYMBOL	KIF26A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcgagagtgggctgcctcccagggcgccgcaccgcagcctcaagtcccccagaagagggccacaggtctgcagcggcgccgctgattccgccccactgcccacaccactgcctgggcccgaagcccagcctcccgggcagtggtggacctgccccgcc
ctggctggctccctgaaggagccgttcgagatcaaggtgtacgagatcgtacgtggagcgcttcagcggccccccccacccgaggaggccccacccaggtctggcgtgcgtcagtaacaagctgcggctggcggagcgcagcagcagcggctgcgggaggtgca
ggccaagcacaagcactgtgtgaggagctggcggagaccagggcgccgtgatgctggagcctggcgtgctggagcagttttaggtggaccggagctggagcccagtcggccgagtagctggcggccctggagcagccacggcgccctggagcagtcgtgaacc
tgtgcaaggcgacgtcatgatggtcacctcttcgacatcagcgttcagccagtgctgccatccccggggccgagggagtgacgtctga