

MRM1-mitochondrial rRNA methyltransferase 1 homolog (S. cerevisiae) Gene

Reference: PTXBC009416

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	rRNA methyltransferase 1, mitochondrial; 16S rRNA (guanosine(1145)-2'-O)-methyltransferase; 16S rRNA [Gm1145] 2'-O-methyltransferase; CTB-75G16.3; mitochondrial large ribosomal RNA ribose methylase; mitochondrial rRNA methyltransferase 1 homolog; mitochon		

For research use only.

GENE DESCRIPTION	mitochondrial rRNA methyltransferase 1 homolog (S. cerevisiae)
GENE ID	79922
ORIGIN SPECIES	Human
ACCESSIONS	BC009416
VECTOR	pDONR223
NCBI SYMBOL	MRM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacgtgttccactgatgacctcaccggattttacagacaaagcccagcagggtggtctgtggccggcagcgtgggtgcccagcacagagatccccagtcctccagatccccatcatgattgcttgagttcctctgggaacggcctactctctgtgctgggaatga
gggtcagggtctatcccaggaggtgagggcctctgcccagcttctcaccatcctgccccggcgcagctgctcctggaacttgagtccttgaacgtctctgtggctgcaggaaattcttctactccatttgcagccagagggaagggtttccccacagaggggagagaaggcagctct
ccaagaccccaagaacctcagccagggtctgaagggtcagcatggctcagcaccaggggtgtcttcaggccagagaaagagaggcaaaatgaggctga