

GNA12-guanine nucleotide binding protein (G protein) alpha 12 Gene

Reference: PTXBC009457

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NNX3; RMP; gep; guanine nucleotide-binding protein subunit alpha-12; WUGSC:H_GS165O14.2; g alpha-12; guanine nucleotide binding protein (G protein) alpha 12; G protein subunit alpha 12		

For research use only.

GENE DESCRIPTION	guanine nucleotide binding protein (G protein) alpha 12
GENE ID	2768
ORIGIN SPECIES	Human
ACCESSIONS	BC009457
VECTOR	pDONR223
NCBI SYMBOL	GNA12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtccgatgaccacgttttaaaggagaaagagagctcctgtggggccctcgggggtggtctcaggtccatttcagctctgcaacagtgacgcgcagcccggtccggagcgtggtgagctttgttgcctctgggtcagcttgcgtgtctctgtgtgttagaatccagagccc
agaggaagtgcagcgggtcctccgccaacggggagagcctcttcgggcgctgttggcgacagcagcgtgtgattcgcgtagcaggggagttgtttgaaacaccttctga