

BARX1-BARX homeobox 1 Gene

Reference: PTXBC009458

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	homeobox protein BarH-like 1; BarH-like homeobox 1; BARX homeobox 1		

For research use only.

GENE DESCRIPTION	BARX homeobox 1
GENE ID	56033
ORIGIN SPECIES	Human
ACCESSIONS	BC009458
VECTOR	pDONR223
NCBI SYMBOL	BARX1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcctggagaaacgcttcgagaagcagaagcttccacgccggacagaatagcttctgagtcctgggcctgagccagttgcaggtgaagacgtggtaccagaatcggaggatgaagtggagaaaatagtgctgcagggcgccggcctggagtcctccaccaagccc
aagggcgccccaagaagaactcaattccaacgagcgagcgagcttactgagcaggagcgcgccaaggtgagagaaaccggcgaggtgccggcgagcccgacagagagccgcgagggactga