

DSTN-destrin (actin depolymerizing factor) Gene

Reference: PTXBC009477

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACTDP; ADF; HEL32; bA462D18.2; destrin; actin-depolymerizing factor; bA462D18.2 (destrin (actin depolymerizing factor ADF) (ACTDP)); epididymis luminal protein 32; destrin, actin depolymerizing factor		

For research use only.

GENE DESCRIPTION	destrin (actin depolymerizing factor)
GENE ID	11034
ORIGIN SPECIES	Human
ACCESSIONS	BC009477
VECTOR	pDONR223
NCBI SYMBOL	DSTN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctcaggagtgcaagtagctgatgaagtatgctgcattttttatgacatgaaagtcgtaaatgctccacaccagaagaatacaagaaaagaaagagctgtcattttttgtctcagtcagacaaaaagtcacattgtagaagaaggcaaaagagatcttggtggagatgttggtga
accataactgatcctttcaagcattttgtgggaatgcttctgaaaaagattgctgctatgctttgatgatgcaagctttgaaacaaaagaatccagaaaagaagattgatgtttttttgtggcaccagaactagcacctctgaaaaagtaaaatgatctatgaagctccaaggatgcaattaaa
aagaaattcaaggcataaaacatgaatgcaagcaatggaccagaagatcgaatcggtgtgtattgctgaaaagttaggtggaatccttaattgtagctttgaaggatgccctgtgtag