



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C4orf16-chromosome 4 open reading frame 16 Gene

Reference: PTXBC009485

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C4orf16; 2C18; GBAR; PRO0971; gamma-BAR; AP-1 complex-associated regulatory protein; adapter-related protein complex 1-associated regulatory protein; gadkin; gamma-1-adaptin brefeldin A resistance protein; gamma-A1-adaptin and kinesin interactor; gammal1-a		

For research use only.

GENE DESCRIPTION	chromosome 4 open reading frame 16
GENE ID	55435
ORIGIN SPECIES	Human
ACCESSIONS	BC009485
VECTOR	pDONR223
NCBI SYMBOL	C4orf16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggaaactgctgtcgtgacgacgtgcttgcgactcttcgcaaggaagcggggcgctgcagcgaatagcggcgcgaggatccaagtatttgaacaactgctcaagaggtgagcacttaacaatagagtttgagaactctagtagaaggtgatgaaggggagagcccggaagc
agtcatagggccttctactgaggagaagaaattgttgacctaaagagaagcgcttatgttccattgccgaaaaacaaaaagatctgataaagaaattcaaaaagagcaagaagggcagagaattgtgcagcaatcatcttccaacatgagagaatacaaaattcaggccagagaatgatc
ttcgaactctgtttgagaataatgaagtcacagtataagttttcgaagtatgactctcaccagatgctacagtgttgacaccaaatcagaagaagcagttgtgatttaagcagcaaaataatacaactatgtggaatgacgacagacatcccttagatctagatggtgggaagtatgaagaagggaa
tgaataaagtcttccaatgagagaagcgttccaaaacagagaagacatctacggcgacacttaagtataccaacaggagactgggaagtaatctgacatcagcctctgatcttccaatgggtcggagtgagtgatattgttgatgccgaatggatgataatggaattccgagt
attctggaattgaactctgtattataaactgttactctgcataaagctgacacagatcagacagatcagacagatcgaagat

More info: Please contact us for any question: contact@proteogenix.fr