

STXBP6-syntaxin binding protein 6 (amisyn) Gene

Reference: PTXBC009499

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC156; amisyn; syntaxin-binding protein 6; syntaxin binding protein 6		
For research use only.			

GENE DESCRIPTION	syntaxin binding protein 6 (amisyn)
GENE ID	29091
ORIGIN SPECIES	Human
ACCESSIONS	BC009499
VECTOR	pDONR223
NCBI SYMBOL	STXBP6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtgcacaaatctctatcagcaaggaaattttgcaccttctgatgaaaggatctgggagctgtccaagtcgaaggaggacaaagaaaagattcctttcttgccaactggaggtcaagggaatatttaacttatctgctgtcagtgacaaacagaaacccacacagggctcca
tcacaaagggtcaaacagtttgaaggctccacatcatttctggagatcacagtggtgctcagcagcttcgccagggttaattggtatcgtatcctaatgggattcggcagagtttgattgtttgaaaatgcttttgaccagtggtgagccagcacagcgtcagaaaaatgcaccttctcca
gatcctccacatacctgccagaggtacacaggaagccagagtttatttaactgccaatccaaaattatgggaggaacagcattcctcattcagctgctgacagcgtgaccagcgcagtgacagaaggcaagccagggccttgaatgagcgtggagagcgattaggccgagca
gaggagaagacagaagacctgaagaacagcgccagcagtttgcagaactgcgcacaagcttgccatgaagcacaatgttga