

RNF185-ring finger protein 185 Gene

Reference: PTXBC009504

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase RNF185; BSK65-MONO1; BSK65-MONO2; BSK65-PANC1; BSK65-PANC2; BSK65-TEST1; BSK65-TEST2; BSK65-TEST3; ring finger protein 185		

For research use only.

GENE DESCRIPTION	ring finger protein 185
GENE ID	91445
ORIGIN SPECIES	Human
ACCESSIONS	BC009504
VECTOR	pDONR223
NCBI SYMBOL	RNF185
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaagcaagggccctcggcctctgcatctcctgagaactccagtgccagggggccagtgaggagcagcaatggcgtgagagcggagggcaggacagcacttgcagtgcaacatctgctggacacagccaagatgccgtcatcagcctgtgtggccacctctctgtt
ggccgtgtttacatcagtggttggagaccagacctaacagacaggtgtgtctgtttgcaaagctggcaltcagccgagacaaggtcatccccctctatggagggggcagcactgggcaacaggaccacagagagaagaccctctctgtctcaaggacagagccagagccggaga
atagagggggatttcaaggatttggattggagatgggtggtccagatgtctttggaattggggcatttcccttgggatatttgcacagcatttaataataatgatggcgccctcctccagctgtccctgggacacccagatgtggacgagcagttcctgtcacgcctctctatttggg
ccctgtgatcatgttctggtcctgattgcctaa