

NCRNA00152-non-protein coding RNA 152 Gene

Reference: PTXBC009508

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NCRNA00152; C2orf59; long intergenic non-protein coding RNA 152		

For research use only.

GENE DESCRIPTION	non-protein coding RNA 152
GENE ID	112597
ORIGIN SPECIES	Human
ACCESSIONS	BC009508
VECTOR	pDONR223
NCBI SYMBOL	NCRNA00152
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacagacacgcgaaatcacgactcagcccccagcagcactctacctgttgcgcgcgatcacagccgggaatgcagctgaaagattccctggggcctgggtccaacgcgccactgtggactctgaggcctctgcatttgcgggtggtctgcctgtga