

HMGN3-high mobility group nucleosomal binding domain 3 Gene

Reference: PTXBC009529

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PNAS-24; PNAS-25; TRIP7; high mobility group nucleosome-binding domain-containing protein 3; TR-interacting protein 7; thyroid hormone receptor interacting protein 7; high mobility group nucleosomal binding domain 3		

For research use only.

GENE DESCRIPTION	high mobility group nucleosomal binding domain 3
GENE ID	9324
ORIGIN SPECIES	Human
ACCESSIONS	BC009529
VECTOR	pDONR223
NCBI SYMBOL	HMGN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgaagagaaagtctccagagaatacagaggcacaagatggatccaaagtaactaaacaggagccacaagacgggtctgccagattgtcagcgaacctgtccacacaaacctgaacccaacaaagaaaaacatctgctaagaagaacacctggagcaagattagcagaggtgctaaagggaagaggaggaagcaggaagctggaaagggaaggcacagaaaactga