

C10orf104-chromosome 10 open reading frame 104 Gene

Reference: PTXBC009530

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C10orf104; APC16; CENP-27; MSAG; bA570G20.3; anaphase-promoting complex subunit 16; centromere protein 27; cyclosome subunit 16; metabolic syndrome-associated protein; anaphase promoting complex subunit 16		

For research use only.

GENE DESCRIPTION	chromosome 10 open reading frame 104
GENE ID	119504
ORIGIN SPECIES	Human
ACCESSIONS	BC009530
VECTOR	pDONR223
NCBI SYMBOL	C10orf104
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgttcacatcatcctcctcagctgggtgggtcagtggaagttctgtcactggatctggttcagtgctcagacctgccccaccacggaagccctttcacctaccccaaggagctggagagatgtagaagatggctctgagagattcctctcgcaatctgttttagctatcaagtggcatccacgttaaacagggtgaacatgatcagcaagttgctcggatggaaaaactagctggttggtagaagagctggaggctgacgagtgcggtttaagcccatcgagcagctgctgggattcacccctctcagggttga