

STRA13-stimulated by retinoic acid 13 homolog (mouse) Gene

Reference: PTXBC009571

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	STRA13; CENP-X; FAAP10; MHF2; centromere protein X; FANCM associated histone fold protein 2; FANCM-interacting histone fold protein 2; Fanconi anemia-associated polypeptide of 10 kDa; retinoic acid-inducible gene D9 protein homolog; stimulated by retinoic		

For research use only.

GENE DESCRIPTION	stimulated by retinoic acid 13 homolog (mouse)
GENE ID	201254
ORIGIN SPECIES	Human
ACCESSIONS	BC009571
VECTOR	pDONR223
NCBI SYMBOL	STRA13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggagcaggagctggatccggtccggaaggagctggtgagcaggctgctgcactgcactcaaggatgacaagacaaagaagcagcagtcctgtggcgtgcggcaggccagcagaagacgcgtccgtgcggacgtggaccagctggagaaggtgcttcgcagctgctcctggactctag