

HIST1H2BC-histone cluster 1, H2bc Gene

Reference: PTXBC009612

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2B.1; H2B/l; H2BFL; dJ221C16.3; histone H2B type 1-C/E/F/G/I; H2B histone family, member L; histone 1, H2bc; histone H2B.1 A; histone H2B.l; histone cluster 1, H2bc; histone cluster 1 H2B family member c		

For research use only.

GENE DESCRIPTION	histone cluster 1, H2bc
GENE ID	8347
ORIGIN SPECIES	Human
ACCESSIONS	BC009612
VECTOR	pDONR223
NCBI SYMBOL	HIST1H2BC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctgagccagccaagtctgtcccgccccgaagaagggtccaagaaggcagtgaccaaaagcgagaaagatggcaagaagcgcaagcgagcgcaagggagagttactctgtgtacgtgtacaagggtgtaaacaggtccatccccgacactggcatcttccaaggc
catgggcatcatgaattcttctgtaacacataatttgagcgcacgcggggcgaggcttccgcctggcgcatatacaacaagcgtcgacccacacccagggagatccagacggccgtgcgctgtcttccggagagctggccaagcacgccgtgtcgagggcaccaaggccg
tcaccaagtacaccagctccaagtaa