

SPAG16-sperm associated antigen 16 Gene

Reference: PTXBC009614

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WDR29; sperm-associated antigen 16 protein; WD repeat domain 29; pf20 protein homolog; sperm-associated WD repeat protein; sperm associated antigen 16		

For research use only.

GENE DESCRIPTION	sperm associated antigen 16
GENE ID	79582
ORIGIN SPECIES	Human
ACCESSIONS	BC009614
VECTOR	pDONR223
NCBI SYMBOL	SPAG16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgtcagcagggatgccagctccggtgagggctctggaagaggcgttggcgtggtttgacggcagccggggacgcgagggacacggcgacgggtggcggctgagggcgctactacctggaacaggtcaccataactgaagcatctgaagatgactatgaat
atgaagagataccagatgacaatttttagcatccagaaggtgaagaagatctggcaaaagcaatcagatggccaagaacaggtctacagatactgaaatttggacggaaaacagttctcttcaagcatgcagctacctgaagtaataagaagactttctgcaatttctgatcaaaatg
ggaatgaccagaactctgattgcttctcagctggaatggtatgagttatatacagaaggagtgactgaacttagaactgttgggaatgttcagatgtctacaccagattatgttttgaaaatgagaacaaaatttaagaagaattgaagcactacaacaagcagctgagtatgtatfttt
taa