

DPH5-DPH5 homolog (S. cerevisiae) Gene

Reference: PTXBC009620

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DPH5 homolog; AD-018; CGI-30; HSPC143; NPD015; diphthine methyl ester synthase; diphthamide biosynthesis methyltransferase; diphthine synthase; protein x 0011; diphthamide biosynthesis 5		

For research use only.

GENE DESCRIPTION	DPH5 homolog (S. cerevisiae)
GENE ID	51611
ORIGIN SPECIES	Human
ACCESSIONS	BC009620
VECTOR	pDONR223
NCBI SYMBOL	DPH5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgcactgtggactgggagaaccattgcattccttgatcatcacaggaggcagcatcacatccaatggagatggagatgctaagtctgtttccataccagaaaatagctcagaatctcaagcatcaatggactttga