



Reference: PTXBC009647

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ALP; PALP; PLAP; PLAP-1; alkaline phosphatase, placental type; alkaline phosphatase Regan isozyme; alkaline phosphomonoesterase; glycerophosphatase; placental alkaline phosphatase 1; alkaline phosphatase, placental		

For research use only.

GENE DESCRIPTION	alkaline phosphatase, placental (Regan isozyme)
GENE ID	250
ORIGIN SPECIES	Human
ACCESSIONS	BC009647
VECTOR	pDONR223
NCBI SYMBOL	ALPP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]