

IL24-interleukin 24 Gene

Reference: PTXBC009681

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C49A; FISP; IL10B; MDA7; MOB5; ST16; interleukin-24; IL-4-induced secreted protein; melanocyte-associated Mda-7; melanoma differentiation-associated gene 7 protein; suppression of tumorigenicity 16 (melanoma differentiation); interleukin 24		

For research use only.

GENE DESCRIPTION	interleukin 24
GENE ID	11009
ORIGIN SPECIES	Human
ACCESSIONS	BC009681
VECTOR	pDONR223
NCBI SYMBOL	IL24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatttcaacagagctgcaaaagcctgtggacttttagccagcagaccctctgcccctcttctgctggcgacagcctctcaaatgcagatgggtgtctccttgcctgggtttaccctgctctctggagccaggtatcaggggcccaggcccaagaattccactttggccctgccaagt
gaaggggtgttccccagaaactgtgggaagccttctgggctgtgaaagacactatgcaagctcaggataacatcacgagtgcccggtgtgctgcagcagaggttctgcagaacgtctcggatgctgagagctgttacctgtccacaccctgctggagttctactgaaaactgtttcaa
aaactaccacaatagaacagttgaagtcaggactctgaagtcattctctactctggccaacaactttgtctcatctgttcacaactgcaaccagtcagaaaaatgagatgtttccatcagagacagtgccacacaggcgtttctgctattccggagacattcaaacagttggacgtagaag
cagctctgaccaaaagcccttgggaagtggaattctctgacctggatgcagaaattctacaagctctga