

NDUFB1-NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 1, 7kDa Gene

Reference: PTXBC009691

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 1, 7kDa
GENE ID	4707
ORIGIN SPECIES	Human
ACCESSIONS	BC009691
VECTOR	pDONR223
NCBI SYMBOL	NDUFB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgaacttactcagattgtcgggaccactgggttcattgttctgtccctatgggattgtcattggatgtatttagacagaaagatgatgaacggcctaactgccttcggacaagaagatgtattttaaaaggaattgcaaccagtgagaagttacctggaagtaa